


## Probe Preparation for *in-situ* Hybridization



### *Design of primers for amplification of DNA fragment of interest.*

1. Copy the sequence of interest and paste it into the box on the webpage: <http://www.bioinformatics.nl/cgi-bin/primer3plus/primer3plus.cgi>.
2. Under the 'General Setting' tab change the values:  
Primer Tm: - Min: 55; Max: 65  
Primer GC%: - Min: 45; Opt: 60; Max: 60  
Max Tm Difference: - 3  
Under 'Advanced Settings' tab  
GC Clamp: 2
3. Click on the 'Pick Primers' tab in the top right corner.
4. Select the best set from the list of primers.

### *Amplification of DNA fragment of interest*

5. Resuspend the lyophilized primers using molecular grade water to make a stock solution of 100 ng/μl. Prepare a working solution of 10 mM as described above.
6. Add the following reagents in a 200 μl PCR tube:

**Table.** Reaction mixture for amplification of DNA of interest.

Reagents	Volume (μl)
2X PCRBIO Taq Mix Red	12.5
Gene F	1
Gene R	1
Template (gDNA/plasmid/cDNA)	1
Molecular grade water	9.5

**▲CRITICAL STEP:** Prepare at least 5 tubes in order to identify the most optimal annealing temperature in a gradient PCR reaction.

7. Setup the gradient PCR reaction with following conditions:

**Table.** PCR condition for amplifying DNA of interest.

Temperature (°C)	Time (secs)	Number of cycles
95	60	1
95	15	40
Gradient (55-65)	15	
72	15	
4	∞	1

8. Run the reaction mixture in 1% agarose gel for 30 mins.

**⏸PAUSE STEP:** The PCR reaction mixture can be stored at 4°C overnight.

### *Cloning of amplified DNA fragments (using pGEM®-T Vector System)*

1. Add the following reagents in a 1.5 ml microcentrifuge tube (ligation mixture):

**Table.** Reaction mixture for ligating DNA of interest to plasmid.

Reagents	Volume (μl)
2X Rapid ligation buffer	5
pGEM®-T vector	0.5

Amplified DNA (200ng/μl)	0.5
T4 DNA ligase	1
Molecular grade water	3

- Incubate the reaction mixture at 4°C for 16 hrs.
- Take out 1 vial of competent cells and keep the tube on ice for 15 mins.
- Transfer 5 μl of ligation mixture into the competent cell tube and tap gently to mix the solution.
- Leave the mixture on ice for 30 mins.
- Heat shock the cells by transferring the tube into a water bath at 42°C for 45 secs.
- ▲ CRITICAL STEP:** Be careful not to exceed the heat shock step above 45 secs.
- Transfer the tube into ice and leave it for 2 mins.
- Add 500 μl of autoclaved LB broth and incubate the cells in bacterial incubation chamber at 37°C with shaking speed of 225 rpm for 2 hrs.
- Centrifuge the tube at 3000 rpm for 4 mins.
- Inside a biological safety cabinet add the following reagents to an LB agar plate:

**Table.** Reagents for the screening of positive bacterial colonies.

Reagents	Volume (μl)
IPTG	25
X-GAL	25
Ampicillin	25

- Spread the reagents on the plate using glass beads and let the plate dry inside the hood.
- Add 50 μl of supernatant from step 9 and spread across the plate using the glass beads.
- Once dried, seal the plate using parafilm and incubate the plate inside a bacterial incubator at 37°C for 14 hrs.

### *Colony PCR on transformed clones*

- In a 1.5 ml microcentrifuge tube add 10 μl molecular grade water. Pick the transformed white colonies and transfer them into the tube. Vortex gently to homogenize the colony.
- ▲ CRITICAL STEP:** Prepare at least 10 clones for testing.
- Add the following reagents in 200 μl PCR tubes:

**Table 21.** Reaction mixture to verify DNA inserted into plasmid.

Reagents	Volume (μl)
2X PCR BIO Taq Mix Red	12.5
M13F primer	1
M13R primer	1
Homogenized clone	1
Molecular grade water	9.5

- Setup the PCR reaction with the following conditions:

**Table.** PCR settings to amplify inserted DNA along with the M13 region of the plasmid.

Temperature (°C)	Time (secs)	Number of cycles
95	60	1
95	15	30
57	15	
72	15	
4	∞	1

- Run the reaction mixture in a 1% agarose gel for 30 mins and note down the colonies with a single band of the expected size, e.g., those that don't have an empty plasmid.
- Inside a laminar hood add 5  $\mu$ l of Ampicillin stock solution into a test tube with 5 ml LB broth. Transfer 5  $\mu$ ls of homogenized cells from step 1. Do this step for every positive colony.
- Incubate the tubes in a bacterial incubation chamber at 37°C and 225 rpm for 14-16 hrs.

***Isolation of plasmids from transformed clones (using GeneJET Plasmid Miniprep Kit)***

- Harvest the cells in a 1.5 ml centrifuge tube at 3000 rpm for 5 mins (pellet can be stored in 40% glycerol at -80°C for future use).
- Discard the supernatant and resuspend the pellet in 250  $\mu$ l of resuspension buffer.
- Add 250  $\mu$ l of lysis buffer and mix by inverting the tube 6-10 times.
- Add 350  $\mu$ l of neutralization buffer and mix by inverting the tube 6-10 times.
- Centrifuge at 14000 rpm for 5 mins and transfer the supernatant to GeneJET spin column.
- Centrifuge the column at 14000 rpm for 30 secs.
- Add 500  $\mu$ l of wash buffer and centrifuge at 14000 rpm for 30 secs. Discard the flow through and repeat this step one more time.
- Centrifuge the empty column at 14000 rpm for 1 min.
- Transfer the column to a 1.5 ml microcentrifuge tube and add 20  $\mu$ l of elution buffer or molecular grade water. Incubate the mixture at room temperature for 3 mins.
- Centrifuge the column at 14000 rpm for 1 min and measure the concentration of plasmid using Nanodrop.
  - PAUSE STEP:** Prepare a working concentration of 100ng/ $\mu$ l. The purified plasmid can be stored at 4°C for over one month. For long-term storage use a -20°C freezer.

***Perform orientation check on the isolated plasmid***

- Perform PCR with following set of primers to see which one is giving band.

1	2	3	4
M13F	M13F	M13R	M13R
Gene F	Gene R	Gene F	Gene R

Orientation 1



You will observe band with Primer set 2 and 3. Use SP6 polymerase for AS probe.

## Orientation 2



You will observe band with Primer set 1 and 4. Use T7 polymerase for AS probe.

### *M13 amplicon preparation for probe*

1. Add the following reagents in 200  $\mu$ l PCR tubes:

**Table.** Reaction mixture to verify DNA inserted into plasmid.

Reagents	Volume ( $\mu$ l)
2X PCRBIO Taq Mix Red	12.5
M13F primer	1
M13R primer	1
plasmid	1
Molecular grade water	9.5

2. Setup the PCR reaction with the following conditions:

**Table.** PCR settings to amplify inserted DNA along with the M13 region of the plasmid.

Temperature ( $^{\circ}$ C)	Time (secs)	Number of cycles
95	60	1
95	15	30
57	15	
72	15	
4	$\infty$	1

Run the reaction mixture in a 1% agarose gel for 30 mins.

### *RNA probe preparation*

1. Add the following reagents in 1.ml tube:

**Table.** Reaction mixture for probe synthesis

Reagents	Volume ( $\mu$ l)
10x RNA polymerase buffer	2.5
Dig Mix	2.5
M13 amplified DNA	3
T7/SP6 polymerase	2
Ribolock	0.5
Molecular grade water	14.5

2. Incubate the tube at 37 $^{\circ}$ C for 2 hrs.

*Purification of RNA probe (via ethanol precipitation)*

3. Add 80  $\mu\text{l}$  of molecular grade water to the reaction tube from the previous step to raise the volume to 100  $\mu\text{l}$ .
4. Add 10  $\mu\text{l}$  of 3M NaOAc, 10 $\mu\text{l}$  Linear acrylamide, and 200  $\mu\text{l}$  of 100% ethanol.
5. Vortex the mixture for 10 secs and store at  $-20^{\circ}\text{C}$  for 15-20 mins.
6. Centrifuge the mixture at  $4^{\circ}\text{C}$ , 14000 rpm for 15 mins.
7. Carefully remove the supernatant.  
**▲ CRITICAL STEP:** Be very careful not to disturb the pellet.
8. Dry the sample in a vacuum concentrator and add 20  $\mu\text{l}$  of molecular grade water.
9. Prepare a stock concentration of 600 ng/ $\mu\text{l}$  by adding additional water (after a Nanodrop reading) and store aliquots at  $-20^{\circ}\text{C}$ .  
**Ⓜ PAUSE STEP:** RNA can be stored at  $-20^{\circ}\text{C}$  for over 1 year.